




3D Glass-Ceramic BioArrays

Microwell - 20K

20,000 Microwells - 125 Micron Diameter - 25 Microns Deep



Available As:

- Substrates
- Functionalized Slides
- Multiplexed Arrays
- Integrated Systems
- 76 x 25 x 1mm Slides

Benefits

- Reduce reagent usage
- Reduce evaporation
- Excellent spot morphology
- High probe densities
- Improved CV's
- Improved dynamic range
- Lower assay costs
- Functionalize with your surface treatment chemistries

