





3D Glass-Ceramic BioArrays

# Micropillar 2700 - 400

2700 Micropillars - 400 Micron Diameter - 200 Microns Tall ( $\pm 25\mu\text{m}$ )



- 2700 Analytes on a Slide
- Excellent Spot Morphology
- Great Solid Phase Kinetics
- No Dye Ratios – More Dyes
- High Throughput
- 1-150 Multiplex/Sample
- No Beads, No Flow Cyto
- Reduced Reagent Use
- Isolated Interrogation Plane
- Flat, Smooth Pillar Top
- Optically Clear Pillars

## Available As:

- Substrates
- Functionalized Slides
- Multiplexed Arrays
- Integrated Systems
- 76 x 25 x 1 mm slides

